STERİL VE TEK KULLANIMLIK KAN GAZI ŞİRINGASI
KULLANIM KILAVUZU /
STERILE AND SINGLE USE BLOOD GAS SYRINGE
USER MANUAL

AÇIKLAMA / DESCRIPTION

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Genject Blood Gas Syringe products; It consists of the main components of the piston, body, gasket and there are options with/without needle.

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Blood Gas syringes manufactured by injection method from high molecular weight polypropylene and polyethylene raw material are assembled with needle and isoprene rubber manufactured from SUS304 stainless material in accordance with ISO 5832-1 & ASTM F138 standard and the final product is obtained.

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The Blood Gas Syringe is primarily indicated for the determination of unbound dissolved oxygen and carbon dioxide in blood and for blood acid/base analysis.

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For measurement of blood gases, arterial blood samples, usually extremity taken from the veins.

When preparing the syringe for use, the syringe package is opened from the "open here" section indicated on the label.

To insert the needle, it is removed from its protective sheath and gently pushed through the conical head of the syringe.

It is checked whether the needle is properly connected to the conical head.

To minimize blood loss, capillary blood may sometimes be drawn instead of arterial blood. In order to ensure vasodilation and adequate blood flow, the area from which blood will be drawn should be kept at 45°C for 10 minutes. After wiping the first drop, spontaneous blood should be drawn into capillary tubes containing lyophilized heparin.

Arterial blood should always be drawn with a syringe and a tourniquet should not be used.

After entering the vein, the plunger of the syringe should not be pulled. Because the blood pressure of the artery will allow the blood to fill the injector.

The most appropriate method is to collect arterial blood with sterile syringes containing 2 or 2.5 mL of lyophilized heparin in an anaerobic environment.

When using liquid heparin, the size of the syringe and the amount of heparin and blood in the syringe are important.

After washing the inside of the syringe with the liquid heparin drawn into the syringe, this heparin should be emptied together with the air inside the syringe. In this way, adequate anticoagulation will be achieved, 0.05 mg of heparin per mL of blood (EDTA, citrate, and other anticoagulants change blood pH).

When the heparin in the syringe is not fully discharged or a small amount of blood is drawn into the syringe; Due to the dilution due to liquid heparin, especially the PCO2 values will be affected and the blood pH will be falsely low.

During blood collection, it is necessary to prevent or minimize the contact of blood with air. Even with extreme care, the blood sample in the syringe may come into contact with air remaining in the needle or neck of the syringe. Immediately after the blood is drawn, the air bubble in the syringe should be expelled slowly and the tip of the needle should be covered with a material that will prevent the passage of air, such as gum, glass paste, soap, stopper; The needle should never be bent. Because, compared to arterial blood, PCO2 in dry air is much lower; whereas PO2 (~155 mm Hg) is about 60 mm Hg higher. Therefore, when the blood taken from a person breathing normal room air comes into contact with the air, the CO2 content of the blood and the PCO2 amount decrease; O2 content will increase.

On the other hand, in people whose blood PO2 level is higher than 150 mm Hg due to O2 treatment, etc., the contact of blood with air will lead to a decrease in the O2 amount.

To minimize blood loss, capillary blood may sometimes be drawn instead of arterial blood. In order to ensure vasodilation and adequate blood flow, the area from which blood will be drawn should be kept at 45°C for 10 minutes. After wiping the first drop, spontaneous blood should be drawn into capillary tubes containing lyophilized heparin.

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Storage and conditions of use shall be such as to guarantee that the device is in perfect condition and packaged. No damaged products should be used. The expiry date given on the label guaranteeing the sterility of the device must be checked. It is important to check the integrity of the package.

STERİL VE TEK KULLANIMLIK KAN GAZI ŞİRINGASI
KULLANIM KILAVUZU /
STERILE AND SINGLE USE BLOOD GAS SYRINGE
USER MANUAL

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This device is supplied sterile. Check the package label for more information.

The device has been sterilized by GAMMA ray sterilization in accordance with current standards. The GAMMA label symbol on the label and/or sterilization package will be purple.

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Devices must be stored in a dry place in their original unopened packaging in a range of 25-5°C, 10-85% RH and must not be used after the expiry date.

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It can be used for all age groups. It is designed for use by healthcare professionals and clinical support specialists.

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To request additional information, please refer to the contact information listed in this document.

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