STERİL BESLENME (KATETER) ENJEKTÖRÜ
KULLANIM KILAVUZU /
STERILE FEEDING (CATHETER) SYRINGE
USER MANUAL

AÇIKLAMA / DESCRIPTION

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The Genject Sterile and Single-Use Catheter Syringe product consists of the plunger and barrel as its main components. It is needle-free and has a volume of 50 ml.

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Feeding syringes manufactured by injection method from high molecular weight polypropylene and polyethylene raw material are assembled with isoprene rubber and the final product is obtained.

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Feeding syringes are used to inject the necessary medication or nutrients with the help of a syringe by making a mouth of the abdominal wall of the stomach in order to ensure the feeding of patients who cannot be fed orally.

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WARNINGS

- Do not use products that are torn or whose package has been opened.
- Open the package from the "Open it here" part on the label.
- After removing the syringe from the package, check for cracks that may cause discharge.
- Make sure there is no air in the syringe.
- Do not reuse.
- Do not resterilize.
- Observe the warnings on the label information before using the syringe.
- Use in accordance with the expiry date.
- Do not use the product that has expired.
- Do not use together with different devices other than the intended use and do not use outside the intended purpose.
- This product must only be used by trained personnel.
- Avoid exposure of content to light.
- Ensure that the feeding nutrition does not come into contact with the outer surface of the product body (the surface with the measurement markings).
- Ensure that the conical tip of the syringe is immersed in the feeding nutrition before drawing it into the syringe.
- Store at 5-25 °C.

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Storage and conditions of use shall be such as to guarantee that the device is in perfect condition and packaged. No damaged products should be used. The expiry date given on the label guaranteeing the sterility of the device must be checked. It is important to check the integrity of the package.

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This device is supplied sterile. Check the package label to learn more. The device has been sterilized by ethylene oxide sterilization (ETO) in accordance with current standards. The ETO label icon on the label and/or sterilization packaging shall be in green color.

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Devices must be stored in a dry place in their original unopened packaging in a range of 25-5°C, 10-85% RH and must not be used after the expiry date.

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It can be used for all age groups. It is designed for use by healthcare professionals and clinical support specialists.

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To request additional information, please refer to the contact information listed in this document.

STERİL BESLENME (KATETER) ENJEKTÖRÜ
KULLANIM KILAVUZU /
STERILE FEEDING (CATHETER) SYRINGE
USER MANUAL

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